

Tarso Clubfoot AFO by MARKELL

A semi-rigid AFO for use in the Ponseti Method of Clubfoot Treatment

Treatment of Idiopathic Clubfoot

The Ponseti Method has become the worldwide standard of care for the treatment of idiopathic clubfoot. It consists of a series of long leg casts and a percutaneous tenotomy, followed by an extended period of bracing to maintain correction.

The Bracing Phase

The all important bracing phase of treatment is implemented after the baby comes out of the last cast. During this phase of treatment, it is essential that the feet be maintained in corrected position, to prevent a relapse from occurring. The Ponseti Method of treatment requires that bracing of the feet in corrected position, be continued until 4 to 5 years of age.

THE TARSO CLUBFOOT AFO PROVIDES THE FOLLOWING FEATURES:

- Symmetrical shape for use on either foot.
- Semi-Rigid Construction with soft lambskin lining
- Convenient reversible hook and loop closures at top and bottom.
- Strong secure buckle and strap with dorsal pad to lock the foot down and hold it in place.
- Rigid base plates with threaded receptacles for mounting to the Dobbs Dynamic Clubfoot Bar, Horton Click Bar, or traditional Denis Browne abduction bars.

TARSO CLUBFOOT AFO MODEL # 1646 IS AVAILABLE
IN NUMBER SIZES FROM 0000 UP TO 9
SUGGESTED I CODE: L1930* X 2

*The final and sole responsibility for the correct coding, within established laws, rules, and standards of products, rests upon the party submitting the claim.



MARKELL SHOE CO., INCORPORATED

504 Saw Mill River Road, P.O. Box 246, Main Station, Yonkers, NY 10702 • Tel: 914-963-2258 • FAX: 914-953-9293